

Incubator 8000 IC/SC/NC Technical Service Bulletin # 21

Re:	Product improvement for membrane keypad
Date:	July 31, 2001
Reference Doc:	Following Service Bulletins from Dräger Medizintechnik GmbH: <ul style="list-style-type: none">• 6141.221 Incubator 8000 SC # 26, March 2001• 6141.222 Incubator 8000 IC # 7, March 2001
Reason:	It has been reported that disinfecting and cleaning liquids may cause malfunctions of the keypad if these liquids are used excessively.
Solution:	Seal openings of the membrane keypad with silicone rubber.
Devices affected:	All units Incubator 8000 IC, Incubator 8000 SC, and Incubator 8000 NC.
When:	Next Service Call or Preventive Maintenance.
Additional Info:	See conversion instructions " Sealing of membrane keypad " P/N 2M22267
Cost:	Covered by Customer, post warrenty
Ordering Info:	Silicon rubber Wacker Elastosil E41 P/N 12 02 537
Distribution:	Dräger Service Personnel and Authorized Service Organizations for NCS products.

If you have any questions, please contact Technical Support by phone at 1-800-543-5047 or by fax at 1-215-721-5789.

Dräger Medical, Inc.
Technical Product Manager